



- Of all ankle injuries, 85% are sprains and 85% of all sprains involve lateral ligaments.
- One result of these sprains is the chronic disability that occurs.
- Statistically, up to 44% of patients may have persistent symptoms one year following the injury.

Patent No. 5,620,413

Perform 8™

Lateral Ankle Stabilizer

Neoprene sock

- provides compression for increased blood flow
- gives support to soft tissue
- fits comfortably

Non-elastic rigid strap

- wraps in a figure 8 configuration to provide ankle support
- provides non-restrictive comfort zone over extensor tendons in anterior area of ankle
- adjusts easily for maximum dorsiflexion of the foot

Pads parallel to the Achilles tendon

- prevent posterior pressure on the Achilles tendon
- redistribute strap pressure
- protect Achilles tendon

Hook fasteners with touch down pad

- assure continuous adhesion to the sock
- lock strap at terminal end
- make adjustment of strap quick and easy

Item #	Size	Men	Women
50051	Small	5-7	6-8
50052	Medium	7-9	8-10
50053	Large	9-11	10-12
50054	X-Large	11-13	12-15

Material Content: 58% Nylon, Polyester