



Helps prevent
and relieve
wrist pain



Aluminum stay
helps maintain
neutral wrist
position
during sleep



Reverses to
fit either
hand

DESIGNED BY AN ORTHOPEDIC SPECIALIST

Therapeutic Treatment for Night

WrisTimer PM wrist support provides therapeutic treatment during the night while allowing a comfortable fit to prevent flexion of the wrist. The aluminum stay holds the wrist in a neutral position to reduce inflammation and provide support.



Item numbers:
60500 (Universal)

Key benefits:

Dorsal application allows a comfortable fit and prevents flexion of the wrist

Maintains neutral wrist positioning during sleep

The aluminum stay holds the wrist in a neutral position

One size fits all

Ideal for daily activities:

Repetitive tasks

Post-op rehab

Carpal tunnel syndrome

Improving patient compliance at night

Material Content: Cotton, Nylon, Nickel, Polyurethane